

CONDUCTOMETRIC TITRATION DETERMINATION OF STRENGTH OF AN ACID

**Dr.HARINDRAN SUHANA
ASSISTANT PROFESSOR
DEPARTMENT OF CHEMISTRY
SRM INSTITUTE OF SCIENCE AND TECHNOLOGY
KATTANKULATHUR**

Aim: To determine the strength of a given solution of HCl using the given NaOH solution by conductometric titration.

- Apparatus required
- 1) Burette
- 2) pipette
- 3) Standard measuring flask
- 4) glass rod
- 5) Beaker
- 6) Wash bottle
- 7) funnel
- 8) Conductivity meter
- 9) conductivity cell

Reagents

- 1) sodium hydroxide solution(0.1N)
- 2) distilled water
- 3) HCl solution whose concentration is to be determined


Principle:

- Conductance depends upon the number of ions present in a solution and their ionic mobility.
- For a neutralization reaction between an acid and a base, the addition of the base would lower conductivity of the solution initially.
- This is because the H^+ ions would be replaced by the slow moving Na^+ ions of the base.
- After the equivalence point is reached, further addition of excess alkali introduces fast moving OH^- ions and hence the conductance of the solution increases.
- Therefore, two straight lines with opposite slopes will be obtained when the conductance values are plotted graphically against volume of sodium hydroxide added.
- The point where these two lines intersect is the equivalence point.

Procedure:

- Burette solution: NaOH solution
- Pipette solution: 10ml of the made up HCl solution
- The given HCl solution is made up to 100ml in a standard flask.
- 10ml of the made up HCl solution is pipetted out into a 250ml beaker.
- It is diluted with 90ml distilled water so that the conductivity cell is completely immersed in the solution.
- The solution is stirred well with a glass rod and the initial conductance of the solution is noted.
- 1ml NaOH is added from the burette and conductance is measured after each addition.
- Initially, the conductance decreases and at one point it starts increasing.
- The point at which the conductance starts increasing shows the end point has been reached.
- A plot of Conductance against volume of sodium hydroxide gives the end point.
- The intersection of the two straight lines and the volume corresponding to the intersection point to the x-axis is noted.
- In order to get the accurate end point, the fair titration is performed in the same manner, by adding NaOH in increments of 0.2ml before and after the end point.
- The strength of HCl can be calculated using the formula
- $V_1N_1 = V_2N_2$

Model Graph(Common for both pilot and fair titration)



Conductometric titration of a strong acid (HCl) vs. a strong base (NaOH)

Table 1: Standard sodium hydroxide Vs unknown Hcl(Pilot titration)

[illegible]

Table 1: Standard sodium hydroxide Vs unknown HCl(Fair titration)

[illegible]

- Calculations
- V_1 = Volume of HCl solution (10ml)
- N_1 = Normality of HCl solution
- V_2 = Volume of NaOH (From the graph)
- N_2 = Normality of NaOH (0.1N)
-
- $N_1 = V_2 N_2 / V_1$
- Hence N_1 - the strength of HCl can be calculated
-
-
- Result:
- The strength of the given HCl solution = _____ N