

S. Kunal Keshan
RA2011004010051

ECE – A

**Physics: Electromagnetic
Theory, Quantum
Mechanics, Waves and
Optics- 18PYB101J**

EXPERIMENT- 4

09

21-05-2024

DETERMINATION OF PARAMAGNETIC SUSCEPTIBILITY QUINCKE'S METHOD.

AIM:

To measure the susceptibility of paramagnetic solution by Quincke's tube method.

APPARATUS REQUIRED:

Quincke's tube, travelling microscope, Sample (FeCl_3 solution), electromagnet, power supply, Gauss meter.

PRINCIPLE:

Based on molecular currents to explain Para and diamagnetic properties. magnetic moment to the molecule and such substances are attracted in a magnetic field are called paramagnetic.

The repulsion of diamagnetic is assigned to the induced molecular current and its respective reverse magnetic moment.

The force acting on a substance, of either repulsion or attraction, can be measured with the help of accurate balance in case of solids or with the measurement of rise in level of narrow capillary in case of liquids.

The force depends on the susceptibility χ , of the material i.e., on ratio of intensity of magnetization to magnetizing field I/H . If the force on the substance and field are measured then value of susceptibility can be calculated.

FORMULA:

The susceptibility of the given sample is found by the formula,

$$\chi = \frac{2(P - \sigma)gh}{H^2} \text{ kg m}^{-1} \text{ s}^{-2} \text{ gauss}^{-2}$$

where P is the density of the liquid or solution (kg m^{-3})
 σ is the density of air (kg m^{-3})

g is the acceleration due to gravity (m/s^2)
 h is the height through which the column rises (cm)
 H is the magnetic field at the centre of pole pieces (Gauss)

OBSERVATION:

ρ = density of the liquid or solution = 1480 kg/m^3

σ = density of air = 1 kg/m^3

For Travelling microscope,

Least Count = 0.001 cm

Microscopic reading without field (h_0) = 2.854 cm

$TR = MSR + (VSC \times LC)$

CALCULATION:

$$MSR + (VSC \times LC) = TR \quad \text{and} \quad h \times 10^{-2} \quad \text{and} \quad h/H^2 =$$

$$1. TR = 2.95 + (2 \times 0.001) \\ = 2.952 \text{ cm}$$

$$; h = 0.02954 \text{ m}; \quad h/H^2 = 6.103 \times 10^{-9}$$

$H = 2200$
↓

$$2. TR = 3.00 + (21 \times 0.001) \\ = 3.021 \text{ cm}$$

$$; h = 0.03021 \text{ m}; \quad h/H^2 = 3.719 \times 10^{-9}$$

$H = 2850$
↓

$$3. TR = 3.05 + (14 \times 0.001) \\ = 3.064 \text{ cm}$$

$$; h = 0.03064 \text{ m}; \quad h/H^2 = 3.189 \times 10^{-9}$$

$H = 3100$
↓

$$4. TR = 3.10 + (11 \times 0.001) \\ = 3.111$$

$$; h = 0.03111 \text{ m}; \quad h/H^2 = 2.335 \times 10^{-9}$$

$H = 3650$
↓

$$\text{Mean } h/H^2 = 3.836 \times 10^{-9}$$

RESULT:

The magnetic susceptibility of the given sample = 1.112×10^{-4}
 $\text{kg m}^1 \text{s}^{-2} \text{ gauss}^{-2}$

DETERMINATION OF PARAMAGNETIC SUSCEPTIBILITY BY QUINCKE'S METHOD




Table: To find the rise in the capillary tube of the solution.

S.No	Current (i)	Field (H)	Travelling microscope Reading (h_1, h_2) cm			Difference $h = h_1 - h_2$ $\times 10^{-2}$ m	$\frac{h}{H^2}$ (cm^{-1})
	Ampere	Gauss	MSR (cm)	VSC (div)	TR (cm)		
1.	1	2200	2.45	2	2.952	0.02952	6.103×10^{-9}
2.	2	2850	3.00	21	3.021	0.03021	3.719×10^{-9}
3.	3	3100	3.05	14	3.064	0.03064	3.188×10^{-9}
4.	4	3650	3.10	11	3.111	0.03111	2.335×10^{-9}

$$\text{Mean } \frac{h}{H^2} = 3.836 \times 10^{-9}$$